


DISPOSABLE NEEDLE GUIDES



- ✓ Used with under ultrasound guided procedures for Transrectal prostate biopsies, Transvaginal biopsies and Transvaginal oocyte retrieval.
- ✓ The disposable designs reduce the risks of cross contamination and infection.
- ✓ All guides accept 16G -18G sizes instruments.
- ✓ Compatible with well-known endocavity ultrasound probes in the market. See table below.

GEOTEK
Healthcare Products

İvedik OSB Ağaç-Metal San. Sit. 1436. Sokak No: 12 06370 Ostim - Ankara / TÜRKİYE
Tel: +90 (312) 395 98 52 • Fax: +90 (312) 395 98 53
e-mail: info@geotekmedical.com • www.geotekmedical.com

GEOTEK[®]
Healthcare Products
CE 1984

DISPOSABLE NEEDLE GUIDES

COMPATIBLE USG PROBES		REFERENCE CODE	NEEDLE GUIDE
HITACHI	EUP-V33, EUP-V33W, EUP-CC331, EUP-CC531, EUP-V53W, EUP-V12	GTK01	
BK	8808	GTK02	
	8818, 8808e	GTK13	
TOSHIBA	PVF-620ST, PVL-625RT, PVM-620ST, PVG-600S,	GTK03	
	PVT-661VT, PVM-651VT, PVQ-641V, PVG-630V, PVN-661VT	GTK04	
	PVM-740-RT, PVT-770-RT	GTK07	
	SSA-790A (AplioXG)	GTK08	
	PVT-781VT, PVU-781VT, PVT-781VTE, PVU-781VTE	GTK130	
	PVU-770-ST (11C4)	GTK14	
ESAOTE	EC123, EC1123, SE3123	GTK08	
	E8-5R10	GTK11	
ALOKA	UST-670P, UST-675P, UST-676P	GTK06	
	UST-984-5, UST-969P-5	GTK100	
	UST-9118, UST-9124	GTK101	
HITACHI/ALOKA	C41V1, EUP-V73W	GTK18	
SHIMADZU	UBR10R-65U	GTK05	
PHILIPS	BP10-5ec, C9-4ec	GTK08-R	
	C8-4V, C9-3V, C9-4V (ClearVue650 & 350)	GTK110	
	C10-4ec	GTK16	

COMPATIBLE USG PROBES		REFERENCE CODE	NEEDLE GUIDE
GE	IC5-9RS	GTK08	
	E8CS, IC5-9-D, E72, E8C, IC5-9H, E721, 6.5MHz, MTZ, E8C-RS	GTK09	
	BE9C, BE9CS	GTK10	
	E7C-RC	GTK11-R	
	RIC-9W-RS (Voluson S), RIC6-12-D, RIC5-9-D, RIC5-9H	GTK140	
ALPINION	IC9-RS (Voluson S10)	GTK141	
	E3-10	GTK08	
SIEMENS	EV3-10, EV3-10H	GTK09-RS	
	BE9-4, EC9-4 (Acuson X-150)	GTK08	
	EC7, EC10C5, 6.5EC10, EC9-4 (Sonoline G20)	GTK11	
	EC9-4W (Acuson S3000), MC9-4	GTK12	
	Endo-V II (Sonoline G50)	GTK120	
SAMSUNG MEDISON	Endo V	GTK121	
	EC9-4 (Sonoline G40, Acuson X150, Acuson X300), EV9-4 (Sonoline G40 Acuson X150, Acuson X300)	GTK122	
	EC4-9IS, EC4-9ES, ER4-9ES	GTK11	
MINDRAY	EV4-9/10ED, ER4-9/10ED, EVN4-9	GTK15	
	V10-4B, V10-4, CB10-4, 6CV1s, 65CV1, 65EC10EA	GTK11	

* Please contact us for your specific requirements or if you can not find the probe in the list.